
ML359

产品编号: D58600

CAS: 1069858-99-6

分子式: C₁₉H₁₈N₄O₃

纯度: 99.16%

InChi Key: InChIKey=GKJQVQCIPPIODP-UHFFFAOYSA-N

Smiles: CCN(C(c1ccco1)=O)c1nn(C)c2nc3cc(OC)ccc3cc12

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO:84 mg/mL (200.99 mM);Ethanol:84 mg/mL (200.99 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: ML359 is a protein disulfide isomerase inhibitor.