
C-Met inhibitor D9

产品编号: D58613

CAS: 299405-67-7

分子式: C₁₇H₁₅N₃O₂

纯度: 97%

InChi Key: InChIKey=QXJDCXDETPSOKS-UHFFFAOYSA-N

Smiles: O=C1NC=2C=CC=CC2C1=NNC(=O)CCC=3C=CC=CC3

外观: 固体粉末

作用通路: Tyrosine Kinase/Adaptors

作用靶点: c-Met/HGFR

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: C-Met inhibitor D9 is a c-Met kinase inhibitor.