
Interferon receptor inducer-1

产品编号: D58641

CAS: 2215120-36-6

分子式: C₁₂H₂₂N₄O₂

纯度: 99.73%

InChi Key: InChIKey=MRKKUGDWODATMF-SECBINFHSA-N

Smiles: CCCC[C@H](CO)NC1=NC(=NC=C1OC)NC

外观: 固体粉末

作用通路: Immunology/Inflammation

作用靶点: IFNAR

溶解性: DMSO:55 mg/mL (216.25 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Interferon receptor inducer-1 Used in the treatment of a disorder in which the induction of interferon is involved. is an interferon (IFN) receptor inducer.