
Aplodan

产品编号: D58664

CAS: 6903-79-3

分子式: C₄H₁₂N₃O₄P

纯度: 99.17%

InChi Key: InChIKey=FOIPWTMKYXWFGC-UHFFFAOYSA-N

Smiles: CN(CCOP(O)(O)=O)C(N)=N

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO:Soluble

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Aplodan exhibits antiischemic and antiarrhythmic capacity and protective action on the cell membrane. Aplodan can be used in research on the ischemic heart or ischemic heart disease with persistent ventricular premature beats.