
3-hydroxy-phthalicanhydrid

产品编号: D58691

CAS: 37418-88-5

分子式: C₈H₄O₄

纯度: 98%

InChi Key: InChIKey=CCTOEAMRIIXGDJ-UHFFFAOYSA-N

Smiles: O=C1OC(C2=C1C=CC=C2O)=O

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 55 mg/ml (335.12 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: 3-Hydroxyphthalic anhydride is used to synthesize 3-hydroxyphthalic anhydride modified human serum albumin.