
Toceranib Phosphate

产品编号: D58703

CAS: 874819-74-6

分子式: C₂₂H₂₈FN₄O₆P

纯度: 98%

InChi Key: InChIKey=AOORBROPMMRREB-HBPAQXCTSA-N

Smiles: OP(O)(O)=O.Cc1[nH]c(\C=C2/C(=O)Nc3ccc(F)cc23)c(C)c1C(=O)NCCN1CCCC1

外观: 固体粉末

作用通路: Angiogenesis; Tyrosine Kinase/Adaptors; Proteases/Proteasome

作用靶点: c-Kit; PDGFR; Tyrosinase; VEGFR

溶解性: DMSO:10mM

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Toceranib phosphate, is a selective inhibitor of the tyrosine kinase activity of several members of the split kinase RTK family, including Flk-1/KDR, PDGFR