
S6K1-IN-DG2

产品编号: D58705

CAS: 871340-88-4

分子式: C₁₆H₁₇BrN₆O

纯度: 98%

InChi Key: InChIKey=CQXAPCMYRSTDGK-UHFFFAOYSA-N

Smiles: BrC1=NNC2=NC=NC(=C12)N3CCN(C=4C=CC=CC4OC)CC3

外观: 固体粉末

作用通路: PI3K/Akt/mTOR signaling

作用靶点: mTOR

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: S6K1-IN-DG2(S6K1InhibitorDG2) is a potent p70S6K inhibitor with an IC₅₀ value of less than 100 nM.