
Adenosine amine congener

产品编号: D58726

CAS: 96760-69-9

分子式: C₂₈H₃₂N₈O₆

纯度: 98.29%

InChi Key: InChIKey=JFRJCQJVFMHZOO-UHFFFAOYSA-N

Smiles: NCCNC(=O)Cc1ccc(NC(=O)Cc2ccc(Nc3ncnc4c3ncn4C3OC(CO)C(O)C3O)cc2)cc1

外观: 固体粉末

作用通路: GPCR/G Protein; Neuroscience

作用靶点: Adenosine Receptor

溶解性: DMSO:31.25 mg/mL (54.20 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Adenosine amine congener (ADAC) is an agonist of selective A1 adenosine receptor,.