
KYP-2047

产品编号: D58738

CAS: 796874-99-2

分子式: C₂₀H₂₅N₃O₂

纯度: 99.7%

InChi Key: InChIKey=SPXFAUXQZWJGCJ-ROUUACIJSA-N

Smiles: O=C(CCCc1ccccc1)N1CCC[C@H]1C(=O)N1CCC[C@H]1C#N

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 25 mg/mL (73.65 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: KYP-2047 is a very potent, selective Prolyl oligopeptidase (POP) inhibitor.