
4-phenoxypiperidine

产品编号: D58764

CAS: 3202-33-3

分子式: C₁₁H₁₅NO

纯度: 98%

InChi Key: InChIKey=KBYPITRKIJKGMD-UHFFFAOYSA-N

Smiles: c1ccc(OC2CCNCC2)cc1

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 27.5 mg/mL (155.15 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: 4-phenoxypiperidines with muscle relaxant and anticonvulsant activities.