
Ac-rC Phosphoramidite

产品编号: D58796

CAS: 121058-88-6

分子式: C₄₇H₆₄N₅O₉PSi

纯度: 99.49%

InChi Key: InChIKey=QKWKXYVKGFKODW-UHFFFAOYSA-N

Smiles: CC(C)N(C(C)C)P(OCCC#N)OC(C(COC(c1ccccc1)(c(cc1)ccc1OC)c(cc1)ccc1OC)OC1N(C=CC(NC(C)=O)=N2)C2=O)C1O[Si](C)(C)C(C)(C)C

外观: 固体粉末

作用通路: Cell Cycle/Checkpoint;DNA Damage/DNA Repair

作用靶点: Nucleoside Antimetabolite/Analog

溶解性: DMSO : 95 mg/mL (105.31 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Ac-rC Phosphoramidite can modify phosphodisulfide oligonucleotides.