
GS87

产品编号: D58798

CAS: 919936-70-2

分子式: C₁₆H₁₁N₅O₂S

纯度: 99.47%

InChi Key: InChIKey=HNTUWEZDVRVSFY-UHFFFAOYSA-N

Smiles: C([S]c1nnc(o1)c2ccncc2)c3nc(no3)c4ccccc4

外观: 固体粉末

作用通路: PI3K/Akt/mTOR signaling; Stem Cells

作用靶点: GSK-3

溶解性: DMSO: 25 mg/mL (74.1 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: GS87 is a highly specific inhibitor of GSK3 (glycogen synthase kinase 3) that induces extensive differentiation of AML cells.