
Pan-RAS-IN-1

产品编号: D58801

CAS: 1835283-94-7

分子式: C₃₆H₄₁Cl₂F₃N₆O₂

纯度: 98%

InChi Key: InChIKey=GMDFJPPJWWJNKF-UHFFFAOYSA-N

Smiles: NC1(CCNCC1)C(NCCN1c(ccc(CN2CCN(Cc(c(Cl)ccc3)c3Cl)CC2)c2)c2c(-
c(cc2)ccc2OC(F)(F)F)c1)=O

外观: 固体粉末

作用通路: GPCR/G Protein; MAPK

作用靶点: MAPK; Ras

溶解性: DMSO : 50 mg/mL (69.67 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Pan-RAS-IN-1 is an inhibitor of pan-Ras. It disrupts the interaction of Ras proteins and their effectors.