
Protein kinase G inhibitor-1

产品编号: D58812

CAS: 354544-70-0

分子式: C₁₄H₁₈N₂O₂S

纯度: 98.11%

InChi Key: InChIKey=PGBKHXTZZAKXMK-UHFFFAOYSA-N

Smiles: O=C(N)C1=C(SC2=C1CCC(C)C2)NC(=O)C3CC3

外观: 固体粉末

作用通路: Tyrosine Kinase/Adaptors

作用靶点: PKA

溶解性: DMSO: 30 mg/mL (107.77 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Protein kinase G inhibitor-1 is a potent Protein kinase G inhibitor, IC₅₀= 0.9 μM.