
Reproxalap

产品编号: D58829

CAS: 916056-79-6

分子式: C₁₂H₁₃ClN₂O

纯度: 99.89%

InChi Key: InChIKey=GUHFUVLKYSQIOQ-UHFFFAOYSA-N

Smiles: OC(C)(C)C1=NC2=CC=C(Cl)C=C2C=C1N

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 55 mg/mL (232.36 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Reproxalap covalently binds aldehydes including malondialdehyde and 4-hydroxynonenal.
Reproxalap is a reactive aldehyde species sequestering agent for the treatment of the dry eye.