

---

## PRT-060318

产品编号: D58842

CAS: 1194961-19-7

分子式: C<sub>18</sub>H<sub>24</sub>N<sub>6</sub>O · 3HCl

纯度: 98.82%

InChi Key: InChIKey=NZNTWOVDIXCHHS-LSDHHAIUSA-N

Smiles: CC1=CC(=CC=C1)NC2=NC(=NC=C2C(=O)N)N[C@@H]3CCCC[C@@H]3N

外观: 固体粉末

作用通路: Angiogenesis; Tyrosine Kinase/Adaptors

作用靶点: Syk

溶解性: H<sub>2</sub>O: 5.6 mg/mL (16.45 mM; Need ultrasonic and warming)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: PRT-060318 (PRT318) is a novel selective inhibitor of the Syk tyrosine kinase, as an approach to HIT treatment.