

---

## MLS0315771

产品编号: D58858

CAS: 727664-91-7

分子式: C<sub>15</sub>H<sub>12</sub>FNOS

纯度: 98.92%

InChi Key: InChIKey=GRGDWLDSPJEJU-UHFFFAOYSA-N

Smiles: Cc1ccc(C)c(c1)-n1sc2cc(F)ccc2c1=O

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 11 mg/mL (40.24 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: MLS0315771 is an effective and biologically active competitive phosphomannose isomerase (MPI) inhibitor (IC<sub>50</sub>: ~1 μM and a K<sub>i</sub>: 1.4 μM).