
Pirenperone

产品编号: D58931

CAS: 75444-65-4

分子式: C₂₃H₂₄FN₃O₂

纯度: 99.78%

InChi Key: InChIKey=HXCNRYXBZNHDNE-UHFFFAOYSA-N

Smiles: Cc1c(c(=O)n2ccccc2n1)CCN1CCC(CC1)C(=O)c1ccc(cc1)F

外观: 固体粉末

作用通路: GPCR/G Protein; Neuroscience

作用靶点: 5-HT Receptor

溶解性: DMSO: 3.67 mg/mL (9.32 mM, Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Pirenperone is a putative 5-HT₂ receptor antagonist with the central antiserotonergic and antidopaminergic action.