
Tozadenant

产品编号: D58958

CAS: 870070-55-6

分子式: C₁₉H₂₆N₄O₄S

纯度: 98%

InChi Key: InChIKey=XNBRWUQWSKXMPW-UHFFFAOYSA-N

Smiles: COc1ccc(N2CCOCC2)c2sc(NC(=O)N3CCC(C)(O)CC3)nc12

外观: 固体粉末

作用通路: GPCR/G Protein; Neuroscience

作用靶点: Adenosine Receptor

溶解性: DMSO:50 mg/mL (123.00 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Tozadenant is an adenosine A_{2A} receptor antagonist(Ki of 11.5 nM and 6 nM on human and rhesus,respectively)