
Zofenopril calcium

产品编号: D50613

CAS: 81938-43-4

保存条件: Store at -20°C for 3 year(Powder); In DMSO or others solvent store at 2-4°C for two weeks, at -20°C for six months.

产品介绍: CAS: 81938-43-4

分子式: $(C_{22}H_{22}NO_4S_2)_2 \cdot Ca$

分子量: 897.17

MDL: MFCD08704648

纯度: $\geq 98\%$ (HPLC)

溶解性: 溶于DMSO

作用靶点: Angiotensin-converting Enzyme (ACE); Reactive Oxygen Species;

作用通路: Metabolic Enzyme/Protease; Immunology/Inflammation; NF- κ B;

产品描述: 佐芬普利钙 (Zofenopril Calcium) 是抗氧化剂和血管紧张素转化酶 (ACE) 抑制剂。