
BAY 41-2272

产品编号: D58991

CAS: 256376-24-6

分子式: C₂₀H₁₇FN₆

纯度: 99.86%

InChi Key: InChiKey=ATOAHNRJAXSBOR-UHFFFAOYSA-N

Smiles: NC1=NC(=NC=C1C1CC1)C1=NN(CC2=C(F)C=CC=C2)C2=NC=CC=C12

外观: 固体粉末

作用通路: GPCR/G Protein

作用靶点: Guanylate cyclase

溶解性: DMSO:36 mg/mL (100 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: BAY 41-2272 is a direct and NO-independent soluble guanylate cyclase (sGC) stimulator.