
HPGDS inhibitor 2

产品编号: D58998

CAS: 2101626-26-8

分子式: C₂₀H₂₄F₂N₂O₃

纯度: 99.82%

InChi Key: InChIKey=YRWMWDLEUYLURJ-SHTZXODSSA-N

Smiles: FC(F)OC1=CC=C2C(N=CC(C(N[C@@H]3CC[C@@H](C(C)(C)O)CC3)=O)=C2)=C1

外观: 固体粉末

作用通路: Immunology/Inflammation

作用靶点: PGE Synthase

溶解性: DMSO:130mg/ml(343.54mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: HPGDS inhibitor 2 is a potent, selective hematopoietic prostaglandin D synthase (h-pgds) inhibitor, IC₅₀ = 9.9 nm.