
CCI-007

产品编号: D59048

CAS: 939228-52-1

分子式: C₁₅H₁₆N₂O₅S

纯度: 97.81%

InChi Key: InChIKey=KGSKTLMBOLPOTC-UHFFFAOYSA-N

Smiles: O=C(OCC)COC(C=CC(C=C1SC(N)=NC1=O)=C2)=C2OC

外观: 固体粉末

作用通路: Apoptosis; Chromatin/Epigenetic

作用靶点: Apoptosis; Histone Methyltransferase

溶解性: DMSO: 25 mg/mL (74.33 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: CCI 007 is a selective inhibitor of MLL-r, CALM-AF10 and SET-NUP214 leukemia with IC₅₀ of 3.5 uM against MLL-r leukemia cell line PER-485