
680C91

产品编号: D59079

CAS: 163239-22-3

分子式: C₁₅H₁₁FN₂

纯度: 98.74%

InChi Key: InChIKey=YBSDQTBCNYWBMX-ONEGZZNKSA-N

Smiles: C1=CC(=CN=C1)/C=C/C2=CNC3=C2C=CC(=C3)F

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 50 mg/mL (209.85 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: 680C91 is a potent and selective inhibitor of TDO, and its K_i is 51 nM.