
Telotristat ethyl

产品编号: D59111

CAS: 1033805-22-9

分子式: C₂₇H₂₆ClF₃N₆O₃

纯度: 98.11%

InChi Key: InChIKey=MDSQOJYHHZBZKA-GBXCKJPGSA-N

Smiles: CCOC(=O)[C@@H](N)Cc1ccc(cc1)-c1cc(O[C@H](c2ccc(Cl)cc2-n2ccc(C)n2)C(F)(F)F)nc(N)n1

外观: 固体粉末

作用通路: Metabolism

作用靶点: Hydroxylase

溶解性: DMSO:98.5 mg/mL (171.31mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Telotristat ethyl is a orally-delivered inhibitor of tryptophan hydroxylase that reduces serotonin production.