

---

## SD-169

产品编号: D59113

CAS: 1670-87-7

分子式: C<sub>9</sub>H<sub>8</sub>N<sub>2</sub>O

纯度: 99.78%

InChi Key: InChIKey=GQMYQEAXTITUAE-UHFFFAOYSA-N

Smiles: NC(=O)c1ccc2[nH]ccc2c1

外观: 固体粉末

作用通路: MAPK

作用靶点: p38 MAPK

溶解性: DMSO: 10 mg/mL (62.43 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: SD 169 is a selective and ATP competitive the MAP kinases p38 $\alpha$  and p38 $\beta$  inhibitor.