
Inarigivir soproxil

产品编号: D59116

CAS: 942123-43-5

分子式: C₂₅H₃₄N₇O₁₃PS

纯度: 98.35%

InChi Key: InChIKey=CJCYTUJOSMYXLE-JDLSZIHUSA-N

Smiles: CO[C@@H]1[C@H](OP(=O)(OC[C@H]2O[C@H](C[C@@H]2O)n2cnc3c(N)ncnc23)SCOC(=O)OC(C)C)[C@@H](CO)O[C@H]1n1ccc(=O)[nH]c1=O

外观: 固体粉末

作用通路: Microbiology/virology; Proteases/Proteasome

作用靶点: HCV Protease

溶解性: DMSO:102 mg/mL (144.96 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Inarigivir soproxil is an agonist of innate immunity. It also displays effective antiviral activity against resistant hepatitis C virus (HCV) variants (EC₅₀s of 2.2 and 1.0 μM for HCV 1a/1b in cells of genotype 1 HCV replicon systems).