
N-(3-Chloro-4-morpholinophenyl)-N'-hydroxyformimidamide

产品编号: D59122

CAS: 339071-18-0

分子式: C₁₁H₁₄ClN₃O₂

纯度: 99.48%

InChi Key: InChiKey=ZTXADXBIGLLOFX-UHFFFAOYSA-N

Smiles: ON=CNc(cc1)cc(Cl)c1N1CCOCC1

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO : 240 mg/mL (938.60 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: TS-011 is a selective 20-Hydroxyeicosatetraenoic acid synthesis inhibitor.