
YAP/TAZ inhibitor-2

产品编号: D59145

CAS: 2762617-31-0

分子式: C₁₉H₁₄F₄N₄O

纯度: 98.84%

InChi Key: InChIKey=GTRCXQAYBUYBBE-UHFFFAOYSA-N

Smiles: O=C(C(F)=C)N1CC(N2N=C(C=3C=CC(=CC3)C(F)(F)F)C=4C=CC=NC42)C1

外观: 固体粉末

作用通路: Stem Cells

作用靶点: YAP

溶解性: DMSO: 90mg/mL (230.6mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: YAP/TAZ inhibitor-2 is a potent and orally active TEAD-YAP/TAZ inhibitor with an EC₅₀ value of 3 nM. YAP/TAZ inhibitor-2 shows anti-proliferative activity. YAP/TAZ inhibitor-2 shows antitumor activity.