

SEW2871

产品编号: D59146

CAS: 256414-75-2

分子式: C₂₀H₁₀F₆N₂O₂S

纯度: 99.62%

InChi Key: InChIKey=OYMNPJXKQVTQTR-UHFFFAOYSA-N

Smiles: O=C1NC(Cc2ccccc2)C(F)(F)F

外观： 固体粉末

作用通路: GPCR/G Protein

作用靶点: LPL Receptor; S1P Receptor

溶解性: DMSO:8.8 mg/mL(20 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Sphingosine-1-phosphate receptor 1 agonist SEW2871 prolongs heterotopic heart allograft survival in mice.