
AZD0095

产品编号: D59153

CAS: 2750001-23-9

分子式: C₂₇H₃₂N₈O₂

纯度: 98.26%

InChi Key: InChIKey=PHDAGSZQGNWGFEXMMPIXPASA-N

Smiles: CC=1N2C(N=C(C)C1O[C@H]3CN(CC3)C4=CC=C(C=C4)C5=CC=C(CN6CCOCC6)N=N5)=NC(C)=N2

外观: 固体粉末

作用通路: Membrane transporter/Ion channel

作用靶点: Monocarboxylate Transporter

溶解性: DMSO: 5.01 mg/mL (10 mM), Need ultrasonic

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: AZD0095 is a selective inhibitor of MCT4 with an IC₅₀ of 1.3 nM.