
Ocaperidone

产品编号: D59159

CAS: 129029-23-8

分子式: C₂₄H₂₅N₄O₂

纯度: 99.94%

InChi Key: InChIKey=ZZQNEJILGNNOEP-UHFFFAOYSA-N

Smiles: Cc1cccn2c1nc(c(c2=O)CCN1CCC(CC1)c1noc2c1ccc(c2)F)C

外观: 固体粉末

作用通路: GPCR/G Protein; Neuroscience

作用靶点: 5-HT Receptor; Dopamine Receptor

溶解性: DMSO:10 mM

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Ocaperidona (Ocaperidone), a very high affinity dopamine D2 antagonist.