
RORy inverse agonist 1

产品编号: D59222

CAS: 529500-72-9

分子式: C₂₂H₂₀F₃N₃O₃S

纯度: 97.31%

InChi Key: InChIKey=BELQBTQMMQYFF-UHFFFAOYSA-N

Smiles: O=C(NCC=1C=CN=CC1)CN(C2=CC=CC(=C2)C(F)(F)F)S(=O)(=O)C3=CC=C(C=C3)C

外观: 固体粉末

作用通路: Others

作用靶点: Others

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: RORy inverse agonist 1 has anti-inflammatory activity and can be used to treat rheumatism and psoriasis.