
Aprocitentan

产品编号: D59237

CAS: 1103522-45-7

分子式: C₁₆H₁₄Br₂N₆O₄S

纯度: 99.17%

InChi Key: InChIKey=DKULOVKANLVDEA-UHFFFAOYSA-N

Smiles: NS(=O)(=O)Nc1ncnc(OCCOc2ncc(Br)cn2)c1-c1ccc(Br)cc1

外观: 固体粉末

作用通路: GPCR/G Protein

作用靶点: Endothelin Receptor

溶解性: DMSO:46 mg/mL (84.22 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Aprocitentan is ETA and ETB antagonist .