
SAR125844

产品编号: D59242

CAS: 1116743-46-4

分子式: C₂₅H₂₃FN₈O₂S₂

纯度: 98.39%

InChi Key: InChIKey=ODIUNTQOXR XOIV-UHFFFAOYSA-N

Smiles: Fc1ccc(cc1)-c1ccc2nnc(Sc3ccc4nc(NC(=O)NCCN5CCOCC5)sc4c3)n2n1

外观: 固体粉末

作用通路: Apoptosis; Tyrosine Kinase/Adaptors

作用靶点: Apoptosis; c-Met/HGFR

溶解性: DMSO:45 mg/mL (81.72 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: SAR125844, a potent and highly selective inhibitor of the MET receptor tyrosine kinase (RTK), for intravenous administration. (IC₅₀ value of 4.2 nmol/L)