

---

## Adenosine 5'-monophosphate monohydrate

产品编号: D50619

CAS: 18422-05-4

保存条件: Store at -20°C for 3 year(Powder); In DMSO or others solvent store at 2-4°C for two weeks, at -20°C for six months.

产品介绍: CAS: 18422-05-4

分子式:  $C_{10}H_{14}N_5O_7P \cdot H_2O$

分子量: 365.24

MDL: MFCD00005750

纯度:  $\geq 99\%$  (HPLC)

溶解性: 溶于DMSO

作用靶点: Adenosine Receptor; HSV;

作用通路: GPCR/G Protein; Anti-infection;

产品描述: Adenosine 5'-monophosphate monohydrate 是一个腺苷 A1 受体 (adenosine A1 receptor) 激动剂。Adenosine 5'-monophosphate monohydrate 对 HSV-1 和 HSV-2 具有显著的抗病毒活性。