
DS-1001b

产品编号: D59272

CAS: 1898207-64-1

分子式: C₂₉H₂₉Cl₃FN₃O₄

纯度: 100%

InChi Key: InChIKey=UPPAWQBZQBNIE-USRGLUTNSA-N

Smiles: CC(C)(C)N.Cc1cn(C(=O)c2c(-c3c(Cl)cc(Cl)cc3Cl)noc2C(C)(C)F)c2cccc(/C=C/C(=O)O)c12

外观: 固体粉末

作用通路: Metabolism

作用靶点: Dehydrogenase; Isocitrate Dehydrogenase (IDH)

溶解性: DMSO:14 mg/mL (22.99 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: DS-1001b is an inhibitor of mutant IDH-1 (Isocitrate Dehydrogenase-1)