
PF429242 dihydrochloride

产品编号: D59385

CAS: 2248666-66-0

分子式: C₂₅H₃₇Cl₂N₃O₂

纯度: 98.67%

InChi Key: InChiKey=GSUZWFZKTIOWTI-MQWQBNKOSA-N

Smiles: O=C(C1=CC=C(C=C1)CN(CC)CC)N([C@H]2CNCC2)CCC3=CC=CC=C3OC.[H]Cl.[H]Cl

外观: 固体粉末

作用通路: GPCR/G Protein

作用靶点: S1P Receptor

溶解性: DMSO: 83.3 mg/mL (165.81 mM); H₂O : 50 mg/mL (103.63 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: PF429242 dihydrochloride can reversibly inhibit S1P (IC₅₀ of 175 nM).