
ERG240

产品编号: D59397

CAS: 1415683-79-2

分子式: C₇H₁₁NaO₃

纯度: 99.12%

InChi Key: InChIKey=QEQCTLAIARYMND-UHFFFAOYSA-M

Smiles: CC(C(C)CCC([O-])=O)=O.[Na+]

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO : 90.0 mg/mL (541.7 mM; Need ultrasonic); H₂O : 90.0 mg/mL (541.7 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: ERG240 is an inhibitor of branched-chain amino acid aminotransferase 1 (BCAT1). ERG240 can be used for the research of rheumatoid arthritis, cancer, and bone disease.