
JNJ-38877618

产品编号: D59400

CAS: 943540-74-7

分子式: C₂₀H₁₂F₂N₆

纯度: 99.52%

InChi Key: InChIKey=KOAWAWHSMVKCON-UHFFFAOYSA-N

Smiles: FC(c1nnc(cc2)n1nc2-c1ccncc1)(c1cc2cccnc2cc1)F

外观: 固体粉末

作用通路: Tyrosine Kinase/Adaptors

作用靶点: c-Met/HGFR

溶解性: DMSO:4.5 mg/mL (12.02 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: JNJ-38877618 is an effective and highly selective inhibitor of Met kinase (IC₅₀s: 2 and 3 nM for wild type and mutant Met, respectively).