
Apararenone

产品编号: D50627

CAS: 945966-46-1

保存条件: Store at -20°C for 3 year(Powder); In DMSO or others solvent store at 2-4°C for two weeks, at -20°C for six months.

产品介绍: CAS: 945966-46-1

分子式: $C_{17}H_{17}FN_2O_4S$

分子量: 364.39

纯度: $\geq 98.5\%$ (HPLC)

溶解性: 溶于DMSO

作用靶点: Mineralocorticoid Receptor;

作用通路: Metabolic Enzyme/Protease;

产品描述: Apararenone (MT-3995) 是新型非甾体类盐皮质激素受体 (mineralocorticoid receptor) 拮抗剂, 正在开发用于治疗糖尿病性肾病和非酒精性脂肪性肝炎。