

bs-19962R

[Primary Antibody]

PRKRIP1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Dog, Horse) Predicted MW.: 21 kDa Subcellular Location: Nucleus
Clonality: Polyclonal		
GeneID: 79706	SWISS: Q9H875	
Target: PRKRIP1		
Immunogen: KLH conjugated synthetic peptide derived from human PRKRIP1: 1-100/184.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: Binds double-stranded RNA. Inhibits EIF2AK2 kinase activity		