

---

## TEMED

产品编号: C05-03008

CAS: 110-18-9

保存条件: Store at RT.

产品介绍: TEMED is used with ammonium persulfate to catalyze the polymerization of acrylamide when making polyacrylamide gels, used in gel electrophoresis, for the separation of proteins or nucleic acids. Although the amounts used in this technique may vary from method to method, 0.1–0.2% v/v TEMED is a "traditional" range.

### PRODUCT SPECIFIC PUBLICATIONS

---

[IF=3.537] Zou Y et al. Comprehensive analysis of human IgG Fc N-glycopeptides and construction of a screening model for colorectal cancer. J Proteomics. 2020 Feb 20;213:103616. Other ; . 31846768