
TRx0237 mesylate

产品编号: D59950

CAS: 1236208-20-0

分子式: C₁₈H₂₇N₃O₆S₃

纯度: 98%

InChi Key: InChIKey=SPCMQFLNOVTUBM-UHFFFAOYSA-N

Smiles: CN(C)c1cc2c(cc1)Nc1c(S2)cc(cc1)N(C)C.CS(=O)(=O)O.CS(=O)(=O)O

外观: 固体粉末

作用通路: Cytoskeletal Signaling

作用靶点: Tau aggregation; Beta Amyloid; Microtubule Associated

溶解性: Soluble in DMSO、H₂O

保存条件: -20°C

产品介绍: TRx 0237 (LMTX) mesylate is a second-generation tau protein aggregation inhibitor for the treatment of Alzheimer's disease (AD) and frontotemporal dementia.