
TP-3654

产品编号: D60025

CAS: 1361951-15-6

分子式: C₂₂H₂₅F₃N₄O

纯度: 99%

InChi Key: InChIKey=XRNVABDYQLHODA-JCNLHEQBSA-N

Smiles: CC(C)(O)[C@H]1CC[C@@H](CC1)Nc1ccc2ncc(-c3cccc(c3)C(F)(F)F)n2n1

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: Pim

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: TP-3654 is a second-generation Pim kinase inhibitor with Ki values of 5 and 42 nM for Pim-1 and Pim-3, respectively.