
TRAP-6

产品编号: D60033

CAS: 141136-83-6

分子式: C₃₄H₅₆N₁₀O₉

纯度: 99%

InChi Key: InChIKey=HAGOWCONESKMDW-UHFFFAOYSA-N

Smiles: CC(C)CC(NC(=O)C(Cc1ccccc1)NC(=O)C(N)CO)C(=O)NC(CC(C)C)C(=O)NC(CCCNC(=N)N)C(=O)NC(CC(N)=O)C(=O)O

外观: 固体粉末

作用通路: GPCR/GProtein

作用靶点: PAR1; Protease-activated Receptor

溶解性: H₂O: 25 mg/mL (33.38 mM; Need ultrasonic)

保存条件: -20°C

产品介绍: TRAP-6 (PAR-1 agonist peptide), a peptide fragment, is a selective protease activating receptor 1 (PAR1) agonist. TRAP-6 activates human platelets via the thrombin receptor. TRAP-6 shows no activity at PAR4