
ML-18

产品编号: D60035

CAS: 1422269-30-4

分子式: C₃₂H₃₅N₅O₅

纯度: 98%

InChi Key: InChIKey=JOKVJNCYOSFDGC-LJAQVGFWSA-N

Smiles: COc1ccc(cc1)C1(CCCCC1)CNC(=O)[C@H](Cc1c[nH]c2ccccc12)NC(=O)Nc1ccc(cc1)[N+](=O)[O-]

外观: 固体粉末

作用通路: GPCR/GProtein

作用靶点: BRS-3; Bombesin Receptor

溶解性: Soluble in DMSO

保存条件: -20°C

产品介绍: ML-18 is a non-peptide bombesin receptor subtype-3 (BRS-3) antagonist with an IC₅₀ of 4.8 μM.