
CC-401(hydrochloride)

产品编号: D60037

CAS: 1438391-30-0

分子式: C₂₂H₂₅ClN₆O

纯度: 99%

InChi Key: InChIKey=OIBVXKYKWOUGAO-UHFFFAOYSA-N

Smiles: Cl.C(CN1CCCCC1)Oc1cccc(c1)-c1n[nH]c2ccc(cc12)-c1ncn[nH]1

外观: 固体粉末

作用通路: MAPK

作用靶点: JNK

溶解性: Soluble in DMSO、H₂O

保存条件: -20°C

产品介绍: CC-401 hydrochloride is a potent inhibitor of all three forms of JNK with K_i of 25 to 50 nM.