
MC-Val-Cit-PAB

产品编号: D60055

CAS: 159857-80-4

分子式: C₂₈H₄₀N₆O₇

纯度: 97%

InChi Key: InChIKey=AMKBTTRWLGVRER-OFVILXPXSA-N

Smiles: CC(C)[C@H](NC(=O)CCCCN1C(=O)C=CC1=O)C(=O)N[C@@H](CCCNC(N)=O)C(=O)Nc1ccc(CO)cc1

外观: 固体粉末

作用通路: Neuroscience

作用靶点: Protease Cleavable; Cleavable; Others; AChR

溶解性: DMSO: 190 mg/mL (331.79 mM; Need ultrasonic)

保存条件: 2-8°C

产品介绍: MC-Val-Cit-PAB is a cathepsin cleavable ADC linker that is used for making antibody-drug conjugate