
SHU-9119

产品编号: D60063

CAS: 168482-23-3

分子式: C₅₄H₇₁N₁₅O₉

纯度: 98%

外观: 固体粉末

作用靶点: MC3R; MC4R; MC5R; Others

溶解性: H₂O: 10 mg/mL (9.31 mM; Need ultrasonic)

保存条件: -20°C

产品介绍: SHU 9119 is a potent human melanocortin 3 and 4 receptors (MC3/4R) antagonist and a partial MC5R agonist; with IC₅₀ values of 0.23, 0.06, and 0.09 nM for human MC3R, MC4R and MC5R, respectively.